

Antimicrobial resistant pathogens in space environment: challenges and potential solutions

Elisabeth Grohmann and Daniela Wischer
Department of Microbiology
Berlin University of Applied Sciences, Germany

Relevance of Antimicrobials in Space Research

- **Astronauts susceptible to infection:**


- Closed system
- Human immune response altered
- Space conditions can increase bacterial virulence, antibiotic resistance, secondary metabolite and EPS production, biofilm formation
- **Damage of technical equipment through bacterial colonisation**
- Common disinfectants often unsuitable (flammable, chemical corrosion)



Application of Antimicrobial Surfaces on the ISS


AGXX[®], GOX and AGXX-GOX-Combinations

- Specially structured antimicrobial coatings
- On diverse carriers, e.g. fleece, plastic, metals



AGXX®

- **Contact catalyst**
 - **Ag** and **Ru** clusters form **microelectrodes**
 - **Microbes are oxidised, O₂ is reduced to ROS**
- ROS cause cell damage / cell death
- AGXX inhibits growth/ kills germs



Clauss-Lendzian et al., 2018


Most Important Results



Ankita Vaishampayan, PhD


Bacterial Load on ISS Was Reduced by 80%

Black: *Staphylococcus* spp.; gray: *E. faecalis*; white: *B. cereus*



Sobisch et al., 2019


AGXX Strongly Reduces Biofilm Formation of MRSA



Vaishampayan et al., 2018

AGXX Reduces Expression of Virulence and Biofilm Genes

9



Vaishampayan et al., 2018

SIRIUS: 4 Month Isolation Experiment in Moscow



Daniela Wischer, PhD

SIRIUS-2019: International Isolation Study

- Crew isolated in Experimental Complex (IBMP, Moscow) for 4 months



AGXX and GOX Exposure in SIRIUS Fitness Room






Approach

- Exposure of AGXX and GOX aboard SIRIUS-2019
- Control experiment in non-isolated environment
- Analyses of bacterial communities from the materials after 1, 2 and 4 months
- Cultivation-based analyses and 16S rRNA gene amplicon sequencing
- Medically relevant isolates were characterised:
 - antibiotic resistance profile
 - biofilm formation
 - plasmid content/transfer potential

Growth Inhibition by Antimicrobials



AGXX- fleece

GOX-fleece

Fleece


Cultured Microbial Communities

- Dominance of *Staphylococcus* spp.
 - Several opportunistic pathogens
 - High biofilm formation capacity
 - Largely MDR
 - No community shift over time
or with materials
-
- Diverse environmental bacteria
 - Fewer potential pathogens
 - Less antibiotic-resistant



16S rRNA Gene Amplicons

- Isolates are not the dominant taxa
- Highly similar communities in both environments
- *Firmicutes*, *Proteobacteria* and *Actinobacteria* > 90%
- *Staphylococcus* sequences comprise only 0.1 – 0.4%





Summary SIRIUS 4 Months Study

- Biofilm-forming, human-derived Staphylococci, opportunistic pathogens, dominate SIRIUS-2019 cultured community
- No resistance to last-resort antibiotics
- Few isolates with large plasmids
- Microbial community stable over 4 months
- 16S rRNA gene amplicons suggest bacterial communities are not dominated by *Staphylococci*



Alternative Experiments to SIRIUS 2022 Study (8 Months)

Experimental Procedures (ISO 22196:2011 with Modifications)

Desiccator with antimicrobial materials adjusted to $\geq 90\%$ relative humidity, materials preincubated to enhance wetting

2 x 2 cm test and control materials

Test strains inoculated during exponential growth (10^9 CFU/mL)

Incubation at 35°C at $\geq 90\%$ relative humidity


70 μL samples were taken at $t = 0, 120, 180$ min, and 24 h

Bacterial survival determined by drop plating of dilutions on LB agar

Incubation for 16 h at 37°C

Survival rate determined as CFU/ cm^2 material

Survival of *E.coli* DSM 682 - High Cell Load



5 log reduction in 2 hours

AGXX-antimicrobial activity in simulated dental care units

Polyamide tubes with 0.1% and 1% AGXX



Test organisms: *Pseudomonas aeruginosa*,
Legionella anisa



Dominique Pütz



CONCLUSIONS

AGXX reduces biofilm formation and consequently antibiotic resistance transfer

AGXX has high antibacterial activity against Gram+ and Gram- bacteria

AGXX demonstrated 8 log reduction of *E. coli* in 150 min

Antimicrobial activity of AGXX is faster than that of Ag

The higher the ruthenium concentration in AGXX the higher the log reduction of the pathogen

THANK YOU



Dr. Daniela Wischer

Dr. Ankita Vaishampayan

Dr. Lydia Sobisch

Dominique Pütz

Robert Hadersbeck

Layla Morad

Dr. Carsten Meyer (Largentec)

Dr. Olaf Wagner (FU Berlin)

Berliner Hochschule für Technik
Studiere Zukunft

